



## TRUSTline Rapid Tests for Dengue



### Dengue IgG/IgM Rapid Test

AR0061C | Serum/Plasma/Whole Blood

- ▶ Detects and differentiates dengue IgG and IgM antibodies
- ▶ Highly sensitive to all 4 serotypes and their isolates
- ▶ IgG: Sensitivity: 100% & Specificity: 98.97%; IgM: Sensitivity: 100% & Specificity: 98.94%
- ▶ Whole blood specimen allows testing at the point of care, saving crucial time
- ▶ Passed MOH registration evaluations in Asia and Latin America
- ▶ No cross-reaction to 17 common diseases, such as HIV, HCV, HBV, CHIK, Malaria, Typhoid
- ▶ No interference with common substances



### Dengue Ag (NS1) Rapid Test

AR0063C | Serum/Plasma/Whole Blood

- ▶ Detects NS1 antigen from all 4 serotypes and their isolates
- ▶ Whole blood specimen allows testing at the point of care, saving crucial time
- ▶ Sensitivity 100% & Specificity: 99.6%
- ▶ Passed MOH registration evaluations in Asia and Latin America
- ▶ No cross-reaction to 17 common diseases, such as HIV, HCV, HBV, CHIK, Malaria, Typhoid
- ▶ No interference with common substances



### Duo Dengue Ag (NS1) + IgG/IgM Rapid Test

AR0062C | Serum/Plasma/Whole Blood

- ▶ Early-stage primary infection: Ag positive
- ▶ Late-stage primary infection: IgM, and/or IgG positive
- ▶ Recent infection: IgG positive
- ▶ Acute secondary infection: Ag and IgG positive and/or IgM Positive

## TRUSTwell ELISA Kits for Dengue



### Dengue IgM ELISA Kit

AE0310 | Serum / Plasma

- ▶ IgM capture ELISA to ensure high specificity
- ▶ Detects IgM for all 4 serotypes or isolates
- ▶ Comparable accuracy to the leading IgM ELISA
- ▶ 2 month open vial stability
- ▶ 18 month shelf life at 2-8 °C storage



### Dengue IgG ELISA Kit

AE0311 | Serum / Plasma

- ▶ IgG capture ELISA to ensure high specificity
- ▶ Detects IgG for all 4 serotypes or isolates
- ▶ Comparable accuracy to the leading IgG ELISA
- ▶ 2 month open vial stability
- ▶ 18 month shelf life at 2-8 °C storage



### Dengue Ag ELISA Kit

AE0312 | Serum / Plasma

- ▶ Detects NS1 for all 4 serotypes or isolates
- ▶ Higher sensitivity than the leading Ag ELISA
- ▶ Detects NS1 as low as 0.3 ng/mL, no hook effect up to 0.1 mg/mL
- ▶ 2 month open vial stability
- ▶ 18 month shelf life at 2-8 °C storage